

FCR Mammography system solution

IP Cassette Reader

**FCR PRIMA Tm**



Table-top type FCR image reader for clinics, boasting the fastest in its class processing speed.

IP Cassette Reader

**FCR CAPSULA XLII**



Compact FCR image reader for hospitals which is suitable for mammography and other purposes, offering high throughput.

IP Cassette

**Type CH**



24 × 30 cm  
18 × 24 cm



24 × 30 cm cassette mounted with the cassette adapter

DRY Imager

**DRYPIX Smart**



Mammography Viewer

**Bellus II**



Analog Mammography unit from FUJIFILM.

the world's first supplier of CR system, offers excellent imaging quality and smooth workflow.

**AMULET FELICIA**  
FUJIFILM ANALOG MAMMOGRAPHY

• The appearance and specifications may be subject to change. • For the details of regulatory information and availability in your country, please contact our local representative.  
• Some items are optional, please contact your subsidiary for the details.  
Manufacturer : METALTRONICA SpA, Distributed by FUJIFILM

# Analog Mammography Unit

AMULET FELICIA with user friendly features, a effective step forward digital mammography in combination with FCR (FUJI Computed Radiography)

## ■ Smart $\mu$ Press compression system

Motor-driven adjustable force and manual compression by means of a double rotary control, are equipped for the aim of appropriate compression of the breast with less discomfort to the patient. Based on the specific density of the breast to be examined, the exclusive FTSE function automatically adjusts the appropriate force to be applied.

## ■ C-arm with isocentric movements

FELICIA's C-arm guarantees isocentric movements which allow all the projections to be carried out with less movement of the patient. The isocentric rotation eliminates the need to adjust the height of the C-arm by changing the projection angle.

## ■ Fully Motorized C-arm

FELICIA is optionally supplied with motorized rotation of C-Arm (preselectable and fine adjustment angles)

## ■ Ultra Software

The microprocessor-controlled automatic exposure meter is calibrated with the advanced ULTRA software. Installation is extremely simple and quick.



### Automatic Collimator

FELICIA can be provided with an automatic collimator. In this case the mammography unit is able to detect the table size (18 x 24 or 24 x 30 cm) and to automatically select the proper collimation field.



### Touchscreen Control panel

Two color displays placed on both sides of C Arm, show information regarding compression force, compressed breast thickness, collimation format, magnification factor, projection angle and patient name. Through a touchscreen the operator can select breast laterality and ACR codes.



### Potter-bucky

There are two different tables (18 x 24 or 24 x 30 cm format) with carbon fiber grid, complete of ID labels. Tables are perfectly interchangeable. (Option)



FUJIFILM supports the Pink Ribbon Campaign by creating and selling Pink Ribbon related goods to its affiliated companies, and donating a part of the sales to "JAPAN SOCIETY OF BREAST HEALTH".

AMULET **FELICIA**