


Small Channel Devices & Accessories




 **STERIS**

Endoscopy

GI Procedure Products

STERIS Endoscopy offers a variety of devices suited for cases requiring endoscopes with smaller channels ranging from 1.0mm to 2.0mm. From foreign body retrieval to tissue sampling and endoscope accessories – we have you covered!


Our specialty bite blocks feature a sturdy, yet forgiving material and smooth insertion rim that help protect delicate instruments and provide comfort for patients.



Pediatric bite block

BX00712808

- » latex-free strap
- » without dental retention rim
- » 38FR



Intermediate bite block

BX00712809

- » latex-free strap
- » without dental retention rim
- » 48FR

Endo-Boot™

BX00711621

- » single-use foam bumper cushions your endoscope and its delicate components during transport and storage helping to reduce the chance of damage

ProArmor™ tip protectors

BX00711623 / BX00711624

- » single-use tip protector featuring a flexible, non-porous material and an open design that allows for air flow, reducing the chances of retaining moisture and endoscope damage

Small Wonder® cleaning brush

BX00711616 / BX00711617

- » specialized, miniature cleaning brush for the smallest endoscope channels (1.00mm - 2.25mm)
- » non-slip, braided wire shaft allows for steady grip


Polypectomy and Tissue Acquisition

World renowned for its reliability and consistency, the Carr-Locke injection needle is a device that can be trusted during everyday procedures as well as difficult maneuvers.

Carr-Locke injection needle – mini

BX00711822

- » designed to operate in tortuous conditions, including duodenoscopes, and is ideal for injection therapy across a broad range of smaller channel endoscopes
- » patented design ensures consistent needle projection every time



STERIS Endoscopy's line of forceps devices address a variety of needs for tissue sampling. Our disposable biopsy forceps feature a Teflon sheath for ease of insertion and are compatible with pediatric bronchoscopes and gastroscopes.

Disposable biopsy forceps

BX00711261 / BX00711262

- » produced with sharpened cutting edges
- » built with stainless steel cups and state-of-art technology
- » disposable forceps are available with or without a spike



Foreign Body Management

Foreign body retrieval can be challenging. Our unique line of products specifically addresses these difficult procedures.



Raptor® grasping device - mini


BX00711178

- » combines alligator and rat tooth capabilities into one device
- » smaller hybrid jaw and catheter ideal for retrieval of foreign bodies in tight positions or procedures requiring an ultra-slim endoscope

Roth Net® retriever - mini

BX00711057 / BX00711168

- » specifically designed for retrieval of a variety of foreign bodies or polyps in smaller channel endoscopes
- » unique design maintains shape and opens to full capacity across multiple passes



Product Information

Procedural Infection Prevention						
part number	description	sheath diameter (mm)	channel size (mm)	bristle diameter(mm)	length (cm)	qty/box
BX00711616	Small Wonder cleaning brush	0.5	1.00-2.25	2.25	160	15
BX00711617	Small Wonder cleaning brush	0.5	1.00-2.25	2.25	230	15
part number	description	compatible scope OD (mm)			qty/box	
BX00711621	small endoscope tip protector	3.0 - 8.8			50	
BX00711623	ProArmor endoscope tip protector - sterile	4.9-15			50	
BX00711624	ProArmor endoscope tip protector - non-sterile	4.9-15			50	
part number	description	opening size	dental retention rim	qty/box		
BX00712808	pediatric bite block	12.6mm / 38FR	no	25		
BX00712809	intermediate bite block	16mm / 48FR	no	100		
Polypectomy & Tissue Acquisition						
part number	description	sheath diameter(mm)	length (cm)	qty/box		
BX00711261	disposable biopsy forceps - oval with spike (sterile)	1.8	160	5		
BX00711262	disposable biopsy forceps - oval (sterile)	1.8	160	5		
BX00711822	Carr-Locke injection needle - 5mm projection / 25 gauge (sterile)	1.8	230	5		
Foreign Body Management						
part number	description	sheath diameter (mm)	approx. net size (cm)	length (cm)	qty/box	
BX00711057	Roth Net retriever - mini	1.8	160	2 x 4.5	5	
BX00711168	Roth Net retriever - mini (sterile)	1.8	160	2 x 4.5	5	
BX00711178	Raptor grasping device - mini (sterile)	1.9	200	-	5	

© 2020 STERIS Endoscopy or its affiliates. All rights reserved.
Unless otherwise indicated, all marks denoted with ® or ™ are registered with the
U.S. Patent and Trademark Office, or trademarks owned by STERIS or its affiliates.

Follow us:    



Endoscopy

STERIS
5960 Heisley Road
Mentor, OH 44060-1834 – USA
Customer Service: 800-548-4873
Fax: 440-639-4495

www.steris.com

761051 Rev. C