



enteroscopy and capsule delivery devices

Enteroscopy devices

US Endoscopy offers a full line of single-use, enteroscopy devices.

Roth Net[®] retriever

BX00711152

sheath diameter 2.5mm

length 350cm

approximate net size 3x6cm

Designed for removal of excised polyps and foreign bodies within the convolutions of the small bowel.

engineered with a strong, reinforced net for extra security

transparent net material allows for continued examination with easy visibility



AdvanCE[®] capsule delivery device

BX00711144

sheath diameter 2.5mm

length 180cm

Allows for direct endoscopic placement of video capsules* in patients who are unable to swallow or pass the capsule through the pylorus.

facilitates endoscopic delivery of video capsules into the stomach or duodenum

enables small intestinal surveillance in challenging patients with dysphagia, pill phobia, gastroparesis, anatomical abnormalities, and more



*Compatible with video capsules ranging from 10.5mm - 11.5mm in diameter and 23.5mm - 26.5mm in length

Carr-Locke injection needle

BX00711811 / BX00711823

sheath diameter 2.5mm
length 230cm
needle projection 5mm
needle gauge 25

BX00711812 / BX00711824

sheath diameter 2.5mm
length 230cm
needle projection 5mm
needle gauge 23

BX00711813 (Pentax 90 series)

sheath diameter 2.5mm
length 230cm
needle projection 5mm
needle gauge 25

BX00711814

sheath diameter 2.5mm
length 230cm
needle projection 4mm
needle gauge 23

The Carr-Locke needle is an ideal choice for injection therapy during enteroscopy procedures.

patented design ensures consistent needle projection every time in most tortuous conditions

luer-lock, spring-loaded handle automatically retracts the needle, reducing the risk of inadvertent needle sticks and scope damage

Articulator™ injection needle

BX00711808

sheath diameter 2.5mm
length 350cm
needle projection 5mm
needle gauge 25

The Articulator™ needle features a kink resistant, stainless steel sheath to pass smoothly and perform in fully articulated endoscopes.



Biopsy forceps



standard
BX0071**1291**

sheath diameter 2.3mm
length 350cm
oval with spike

hot
BX0071**1295**

sheath diameter 2.5mm
length 350cm
oval

State-of-the-art manufacturing processes ensure that these biopsy forceps are sharp and well aligned, allowing for precise tissue samples every time.

Short Throw™ polypectomy snare

BX0071**1101**

sheath diameter 2.4mm
length 350cm
approximate snare size 25 x 54mm

The Short Throw™ snare offers superior loop memory and the comfort of one-handed deployment.

short throw handle offers a 25% stroke travel reduction

designed to work in tortuous conditions



Bipolar hemostasis probes

BX00711841

size 7FR

length 350cm

dual plug

BX00711843

size 7FR

length 350cm

single plug

BX00711847

size 10FR

length 350cm

single plug

The hemostasis probe features a unique, Hydromer tip to help prevent tissue from sticking to the probe.

Jumbo channel cleaning brush

BX00711606

sheath diameter 1.7mm

length 350cm

bristle diameter 7-9mm taper

Designed with a textured, non-slip catheter and soft, specially designed nylon bristles for maximum cleaning effort in scopes that range from 2.8mm and larger.

Valve/control head cleaning brush

BX00711610 / BX00711611

bristle head length 4cm

bristle diameter 10mm

A full head of soft bristles to aid in the removal of gross material, and a grooved, barbell design for a comfortable, non-slip grip.

For more devices available for single balloon and double balloon enteroscopes (230cm length), please visit www.usendoscopy.com.





*listening... and delivering solutions**

5976 Heisley Road
Mentor, OH 44060
phone (1) 440 / 639-4494
fax (1) 440 / 639-4495
customer service (1) 800 / 769-8226
email global@usendoscopy.com

www.usendoscopy.com

Follow us:    

760862 Rev. B