

Roth Net[®] retriever - mini polyp




Polyp retrieval has never been so easy...

The Roth Net® retriever – mini polyp features a smaller net design, offering enhanced maneuverability and visualization while retrieving polyps within the small confines of the colon.



The Roth Net® retriever – mini polyp:

- offers a smaller snare size designed for polyp retrieval
- features a round snare shape that improves maneuverability and allows for a clear view during retrieval
- contains memory points in the snare that enable the net to open to full capacity and maintain a consistent shape during multiple uses



Roth Net® retriever - mini polyp
actual size (2.5 x 4.5cm)

For more information on the Roth Net® retriever - mini polyp, contact US Endoscopy today.

product	description	sheath diameter (mm)	catheter length (mm)	net size (cm)	sterile	units/box
BX00711158	Roth Net® retriever - mini polyp	2.8	230	2.5 x 4.5	no	10
BX00711159	Roth Net® retriever - mini polyp	2.8	230	2.5 x 4.5	yes	5

Distributed by:



listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

761015 Rev. B

US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

Follow us:    