

Roth Net[®] retriever - select



Roth Net® retriever – select

As the leader in foreign body and polyp retrieval for over 20 years, US Endoscopy is proud to introduce the newest addition to the extensive Roth Net® retriever portfolio – the Roth Net® retriever – select.

The new Roth Net® retriever – select features:

An alternative weave design, offering enhanced durability and strength

Unique netting for improved visibility throughout the procedure



US Endoscopy and our distributor partners are your trusted source for retrieval situations. To learn more about our in-servicing and continuing education offerings, visit www.usendoscopy.com/education.

For more information on the Roth Net® retriever – select, contact US Endoscopy today.

product number	sterile	sheath diameter (mm)	length (cm)	net size (cm)	units/box
BX00711169	yes	2.5	230	3.0 x 6.0	10
BX00711170	no	2.5	230	3.0 x 6.0	5

Distributed by:



listening...and delivering solutions*

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

©2017 US Endoscopy, All rights reserved.
US Endoscopy is a subsidiary of STERIS Corporation.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

Follow us:    